

Delta Series Adjustable LCD Pole Mount Arm

The Ergomounts EM2260W is a highly adjustable, durable LCD TV / monitor mounting solution that can be mounted on either vertical or horizontal poles with diameters from 35mm to 50mm.

The arm has the capacity to extend up to 400mm and fold to within 120mm of the pole. The double articulating design allows the monitor to turn a full 180° left or right. This touch screen stable solution provides 180° of monitor tilt adjustment and can be tightened off in position for stability. Available in a white finish the mount can support monitors up to 10kg in weight (approx 24" monitor). Cables are routed along each arm section. This is an ideal mounting solution for dental environments as it will adjust to fit the majority of leading manufacturers dental chairs.

Part Number: EM2260W

- Durable, adjustable pole mount arm for LCD monitors and LCD TVs.
- Folds to within 120mm of pole and extends up to 400mm.
- Tilt, swivel and portrait to landscape movement.
- Touch screen stable.
- Provides +/-90° tilt and +/-90° angle left to right adjustment.
- Attaches to horizontal and vertical poles between 35mm and 50mm diameter.
- Supplied with 75mm x 75mm and 100mm x 100mm VESA interface plates.
- Can be tightened off in position.

Specifications

Arm Ext.:	120mm to 400mm
Monitor Height Adjust	Dependent on pole
Tilt Adjustment:	+/-90°
Monitor Swivel	+/-90°
Monitor Orientation	Landscape to Portrait
Colour	White (silver interface)
Touch Screen Stable	Yes
Keyboard and Mouse Tray	No
Max Monitor Weight	10kg
Product Warranty	1 Year

T +44 (0)1252 333326 F +44 (0)1252 330006
E sales@ergomounts.co.uk

www.ergomounts.co.uk www.medicalmounts.co.uk



24"



up to 400mm



Dependent on pole



+/-90°



+/-90°



L/P



White



Yes



No



10kg



1 Year



VESA 75x75mm
100x100mm



see it your way...